

UBE2L6 Rabbit pAb

货号: B19239

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:100 - 1:500 IF/ICC: 1:50 - 1:200
理论分子量	10kDa/17kDa
实测分子量	18KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Rat kidney
细胞定位	cytosol,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 50 to the C-terminus of human UBE2 L6 (NP_004214.1).
序列	AFNLRISFPPEYPFKPPMIKFTTKIYHPNVDENGQICLPIISSENWKPCTKTCQVLEALNVLVNRPNIREPLRMDLADLLTQNPE LFRKNAEEFTLRFGVDRPS

靶点信息

研究背景	The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or s hort-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-ac tivating enzymes (E1s), ubiquitin-conjugating enzymes (E2s) and ubiquitin-protein ligases (E3s). This gen e encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme is highly similar in prim ary structure to the enzyme encoded by the UBE2L3 gene. Two alternatively spliced transcript variants en coding distinct isoforms have been found for this gene.
------	--

基因ID	9246
基因名	UBE2L6
Swiss	O14933
别名	UBE2L6;RIG-B;UBCH8

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn