

SOX14 Rabbit pAb

货号: B19241

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	26kDa
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 91-240 of human SOX14 (NP_004180.1).
序列	RYVFPLPYLGDTDPLKAAGLPVGASDGLLSAPEKARAFPPASAPYSLLDPAQFSSSAIQKMGEVPHTLATGALPYASTLGY QNGAFGSLSCPSQHTHTHPSTNPGYVPCNCTAWSASTLQPPVAYILFPGMTKTGIDPYSSAHATAM

靶点信息

研究背景	This intronless gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins. Mutations in this gene are suggested to be responsible for the limb defects associated with blepharophimosis, ptosis, epicanthus inversus syndrome (BPES) and Mobius syndrome.
------	--

基因ID	8403
基因名	SOX14
Swiss	O95416
别名	SOX14;SOX28

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn