

SPDYA Rabbit pAb

货号: B19264

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	33kDa/36kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-313 of human SPDYA (NP_877433.2).
序列	MRHNQMCCETPPTVTYVYKSGSNRSHQPKKPITLKRICKDNWQAFEKNTHTNNNKSRRPKGPCLVIQRQDMTAFKLFDDLIQDFLWMDCCCADIADKYLLAMTFVYFKRAKFTISEHTRINFFIALYLANTVEEDEEETKYEIFPWALGKNWRKLFNPKLRLDQLWDRIDYRAIVSRRCCEEVMAIAPHYIWQRERSVHHSGAVRNYNRDEVQLPRGPSATPVDCSLCGKKRRYVRLGLSSSSSLSSHTAGVTEKHSQDSYNSLSMDIIGDPSQAYTGSEVVNDHQSNKGKKTNFLKKDKSMEWFTGSEE

靶点信息

研究背景	Regulates the G1/S phase transition of the cell cycle by binding and activating CDK1 and CDK2. Contributes to CDK2 activation without promoting CDK2 phosphorylation, by inducing a conformation change of the CDK2 T-loop that obstructs the substrate-binding cleft prior to kinase activation. Mediates cell survival during the DNA damage process through activation of CDK2.
------	--

基因ID	245711
基因名	SPDYA
Swiss	Q5MJ70
别名	SPDYA;RINGO3;RINGOA;SPDY1;SPY1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn