

# LSM4 Rabbit pAb

货号: B19271

产品信息

|       |                                                                                                                |
|-------|----------------------------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                                          |
| 宿主    | Rabbit                                                                                                         |
| 克隆性   | Polyclonal                                                                                                     |
| 预测反应  |                                                                                                                |
| 应用    | WB IHC IF/ICC IP                                                                                               |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200<br><b>IP:</b> 1:20 - 1:50 |
| 理论分子量 | 15kDa                                                                                                          |
| 实测分子量 | 14kDa                                                                                                          |
| 形式    | Liquid                                                                                                         |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.         |
| 偶联物   | Unconjugated                                                                                                   |
| 阳性对照  | HepG2,MCF-7,Jurkat                                                                                             |
| 细胞定位  | Nucleus                                                                                                        |
| 纯化    | Affinity purification                                                                                          |

抗原信息


|      |                                                                                                                                               |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-139 of human LSM4 ( NP_036453.1).                             |
| 序列   | MLPLSLLKTAQNHMPMLVELKNGETYNGHLVSCDNWMNINLREVICTSRDGDKFWRMPECYIRGSTIKYLRIPDEIIDMVKEE VVAKGRGRGGLQQKQKQKGRGMGGAGRGVFGGRGRGGIPGTGRGQPEKKPGRQAGKQ |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of the LSm family of RNA-binding proteins. LSm proteins form stable hetero mers that bind specifically to the 3'-terminal oligo(U) tract of U6 snRNA and may play a role in pre-mRNA splicing by mediating U4/U6 snRNP formation. Alternatively spliced transcript variants encoding multiple i soforms have been observed for this gene. |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                  |
|-------|------------------|
| 基因ID  | 25804            |
| 基因名   | LSM4             |
| Swiss | Q9Y4Z0           |
| 别名    | LSM4;GRP;YER112W |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)