

RPH3A Rabbit pAb

货号: B19273

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	76kDa
实测分子量	77kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Cell junction,Membrane,Peripheral membrane protein,synapse
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human RPH3A (NP_055769.2).
序列	MTDTVFSNSSNRWMYPSDRPLQSKLQAGWSVHPGGQPDRQRKQEELTDEEKEIINRVIARA EKMEEME QERIGRLVDRL ENMRKNVAGDGVNRCILCGEQLGMLGSACVVCEDCKKNVCTKCGVETNNRLH SVWLCKICIEQREVWKRS GAWFFKGF PKQVLPQPMPIKTKPQQPVSEPA APEQPAPEPKH PARAPARGDSEDRRGPGQKTGPD PASAPGRGNYGPPVRRASEAR MSSSSRDSESWDHSGGAGDSSRSPAGLRRANSVQASRPAPGSVQ

靶点信息

研究背景	The protein encoded by this gene is thought to be an effector for RAB3A, which is a small G protein that acts in the late stages of neurotransmitter exocytosis. The encoded protein may be involved in neurotransmitter release and synaptic vesicle traffic.
------	--

基因ID	22895
基因名	RPH3A
Swiss	Q9Y2J0
别名	RPH3A

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn