

TKT Rabbit pAb

货号: B19285

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC
推荐浓度	IHC: 1:50 - 1:200
理论分子量	67kDa/68kDa
实测分子量	67kDa/68kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat cells, rat liver tissue
细胞定位	cytosol,extracellular exosome,nuclear body,nucleoplasm,peroxisome
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human TKT (NP_001128527.1).
序列	MESYHKPDQQKLQALKDTANRLRISSIQATTAAGSGHPTSCCSAAEIMAVLFFHTMRYKSQDPRNPHNDRFVLSKGHAAP ILYAVWAEAGFLAEALLNRKISSDLDGHPVPKQAFSTDVATGSLGQGLGAACGMAYTGKYFDKASYRVYCLLDGDELSE GSVWEAMAFASIYKLDNLVAILDINRLGQSDPAPLQHQMIDIYQKRCEAFGWHAIIVDGHVSVEELCKAFGQAKHQPTAIIAK TFKGRGITGVEDKESWHGKPLPKNMAEQIIQEISQIQS

靶点信息

研究背景	This gene encodes a thiamine-dependent enzyme which plays a role in the channeling of excess sugar ph osphates to glycolysis in the pentose phosphate pathway. Alternatively spliced transcript variants encodi ng multiple isoforms have been observed for this gene.
基因ID	7086

基因名	TKT
Swiss	P29401
别名	TKT;HEL-S-48;HEL107;SDDHD;TK;TKT1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn