

ST3GAL4 Rabbit pAb

货号: B19287

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	36kDa/37kDa/38kDa
实测分子量	34-38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse ovary,Rat liver
细胞定位	Golgi apparatus,Golgi stack membrane,Secreted,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 27-267 of human ST3GAL4 (NP_006269.1).
序列	REDSFYFPIPEKKEPCLQGEAASKLFGNYSRDQPIFLRLEDYFWVKTPSAYELPYGTKGSEDLLLRVLAITSSSIPKNIQS LRCRRCVVVGNGHRLRNSSLGDAINKYDVVIRLNNAPVAGYEGDVGSKTTMRLFYPESAHFDPKVENNPDTLLVLVAFKA MDFHWIETILSDKKRVRKGFWKQPPLIWDVNPQIRILNPFFMEIAADKLLSLPMQQPRKIKQKPTTGLLAITLALH

靶点信息

研究背景	This gene encodes a member of the glycosyltransferase 29 family, a group of enzymes involved in protein glycosylation. The encoded protein is targeted to Golgi membranes but may be proteolytically processed and secreted. The gene product may also be involved in the increased expression of sialyl Lewis X antigen seen in inflammatory responses. Multiple transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	6484
基因名	ST3GAL4
Swiss	Q11206
别名	ST3GAL4;CGS23;NANTA3;SAT3;SIAT4;SIAT4C;ST-4;ST3GalA.2;ST3GalIV;STZ;gal-NAc6S

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn