

CHD2 Rabbit pAb

货号: B19311

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IF/ICC
推荐浓度	IF/ICC: 1:50 - 1:200
理论分子量	56kDa/200kDa/211kDa
实测分子量	56kDa/200kDa/211kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 242-501 of human CHD 2 (NP_001036037.1).
序列	SDDLIENTGEGVDEQQDNSETIEKVLDSRLGKKGATGASTTVYAIEANGDPSGDFDTEKDEGEIQYLIKWKGWSYIHSTW ESEESLQQQKVKGGLKLENFKKKEDEIKQWLGKVSPEDVEYFNCQQELASELNKQYQIVERVIAVKTSKSTLGQTDPAHS RKPAPSNEPEYLCKWMGLPYSECSWEDEALIGKKFQNCIDSFHSRNNSKIPTRECKALKQRPRFVALKKQPAYLGGENLE LRDYQLEGLNWLHWSWCK

靶点信息

研究背景	The CHD family of proteins is characterized by the presence of chromo (chromatin organization modifier) domains and SNF2-related helicase/ATPase domains. CHD genes alter gene expression possibly by modification of chromatin structure thus altering access of the transcriptional apparatus to its chromosomal DNA template. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
------	--

基因ID	1106
基因名	CHD2
Swiss	O14647
别名	CHD2;EEOC

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn