

IP6K2 Rabbit pAb

货号: B19332

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IF/ICC
推荐浓度	IF/ICC: 1:50 - 1:200
理论分子量	8kDa/9kDa/11kDa/20kDa/21kDa/49kDa
实测分子量	8kDa/9kDa/11kDa/20kDa/21kDa/49kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human IP6K2 (NP_057375.2).
序列	MSPAFRAMDVEPRAKGVLLPEFVHQVGGHSCVLRFNETTLCCKPLVPREHQFYETLPAEMRKFTPQYKGVVSVRFEEDDRNLCLIAYPLKGDHGIVDIVDNSDCEPKSKLLRWTTNKKHHVLETEKTPKDWVRQHRKEEKMKSHKLEEEFEWLKKSEVLYYTVEKKGNISSQLKHYNPWS

靶点信息

研究背景	This gene encodes a protein that belongs to the inositol phosphokinase (IPK) family. This protein is likely responsible for the conversion of inositol hexakisphosphate (InsP6) to diphosphoinositol pentakisphosphate (InsP7/PP-InsP5). It may also convert 1,3,4,5,6-pentakisphosphate (InsP5) to PP-InsP4 and affect the growth suppressive and apoptotic activities of interferon-beta in some ovarian cancers. Alternative splicing results in multiple transcript variants encoding different isoforms.
------	---

基因ID	51447
基因名	IP6K2
Swiss	Q9UHH9
别名	IP6K2;IHPK2;PIUS

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn