

DCAF7 Rabbit pAb

货号: B19346

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	16kDa/38kDa
实测分子量	16kDa/38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-342 of human DCAF7 (NP_005819.3).
序列	MSLHGKRKEIYKYEAPWTVYAMNWSVRPDKRFRLLALGSFVEEYNNKVQLVGLDEESSEFICRNTFDHPYPTTKLMWIPDT KGVYPDLLATSGDYLRVWRVGETETRLECLLNNKNNSDFCAPLTSFDWNEVDPYLLGTSSIDTTCTIWGLETGQVLGRVN LVSGHVKTLIAHDKEVYDIAFSRAGGGRDMFASVGADGSVRMFDLRHLEHSTIIYEDPQHHPLLRLCWNKQDPNYLAT MAMDGMEVVILDVRVPCTPVARLNNHRACVNGIAWAPHSSCHICTAADDHQALIWDIQMPRAIEDPILAYTAEGEINNVO WASTQPDWIAICYNNCLEILRV

靶点信息

研究背景	This gene encodes a protein with multiple WD40 repeats which facilitate protein-protein interactions and thereby enable the assembly of multiprotein complexes. This protein has been shown to function as a scaffold protein for protein complexes involved in kinase signaling. This highly conserved gene is present in eukaryotic plants, fungi, and animals. The ortholog of this gene was first identified in plants as a key regulator of anthocyanin biosynthesis and flower pigmentation. Alternative splicing results in multiple transcript variants.
基因ID	10238
基因名	DCAF7
Swiss	P61962
别名	DCAF7;AN11;HAN11;SWAN-1;WDR68

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn