

TTC5 Rabbit pAb

货号: B19392

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	48kDa
实测分子量	49kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,MCF7,Mouse spleen,Mouse testis,Mouse thymus,Rat brain
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TTC5 (NP_612385.2).
序列	MMADEEEEEVKPILQKLQELVDQLYSFRDCYFETHSVEDAGRKQQDVQKEMEKTQQMEEVGSGVQGKAQVLMGTGKAL NVTPDYSPKAEELLSKAVKLEPELVEAWNQLGEVYWKKGDVAAAHTCFSGALTHCRNKVSLQNLSMVLRLRTDTEDEH SHHVMDSVRQAKLAVQMDVHDGRSWYILGNSYLSLYFSTGQNPKISQQALSAYAQAQEKVDRKASSNPDLHLNRATLHK YEESYGEALEGFSRAAALDPAWPEPRQREQQLLEFLDRLTSLLESKGGKVKTKKLQSMGLSLRPAH

靶点信息

研究背景	Adapter protein involved in p53/TP53 response that acts by regulating and mediating the assembly of mu lti-protein complexes. Required to facilitate the interaction between JMY and p300/EP300 and increase p5 3/TP53-dependent transcription and apoptosis. Prevents p53/TP53 degradation by MDM2 (By similarity.
基因ID	91875

基因名	TTC5
Swiss	Q8N0Z6
别名	TTC5;Strap

产品验证



Western blot analysis of TTC5 expressed in Mouse spleen,Rat brain,A-431 using TTC5 Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blo cking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn