

LARP1B Rabbit pAb

货号: B19407

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	14kDa/24-38kDa/90-105kDa
实测分子量	105KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,A-431
细胞定位	cytoplasmic stress granule,cytosol,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 170-320 of human LARP 1B (NP_060548.2).
序列	FDYSYGYQEHGERTDQPFQTELNTSMYYDDGTGVQVYPVEEALLKEYIKRQIEYYFSVENLERDFFLRGKMDEQGFLPISLIAGFQRVQALTTNLTNLILEALKDSTEVEIVDEKMRKKIEPEKWPIPGPPPSVPPTDFSQIDCPEFV

靶点信息

研究背景	This gene encodes a protein containing domains found in the La related protein of Drosophila melanogaster. La motif-containing proteins are thought to be RNA-binding proteins, where the La motif and adjacent amino acids fold into an RNA recognition motif. The La motif is also found in proteins unrelated to the La protein. Alternative splicing has been observed at this locus and multiple variants, encoding distinct isoforms, are described. Additional splice variation has been identified but the full-length nature of these transcripts has not been determined.
------	---

基因ID	55132
基因名	LARP1B
Swiss	Q659C4
别名	LARP1B;LARP2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn