

# PIWIL3 Rabbit pAb

货号: B19424

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Mouse,Rat                                                                                            |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                            |
| 理论分子量 | 101kDa                                                                                               |
| 实测分子量 | 101kDa                                                                                               |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Mouse testis,Rat testis                                                                              |
| 细胞定位  | Cytoplasm                                                                                            |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                                                                                                                                                                                                    |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 40-320 of human PIWIL3 (NP_001242904.1).                                                                                                                                                                             |
| 序列   | QLQSTPRPLQEEVPVVRPLQPRAARGGAGGGAQSQGVKEPGPEAGLHTAPLQERRIGGVFQDLVVNTRQDMKHVKDSK<br>TGSEGTVVQLLANHFRVISRPQWVAYKYNVDYKPDIEDGNLRTILLDQHRRKFGERHIFDGNLLLLSRPLKERRVEWLSTT<br>KDKNIVKITVEFSKELTPTSPDCLRYNILFRRTFKLLDFEQVGRNYTKKKAIQLYRHGTSLEIWLGYVTSVLQYENSITLCAD<br>VSHKLLRIETAYDFIKRTSAQAQTGNIREEVTNKLG |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of the PIWI subfamily of Argonaute family proteins. This subfamily of proteins contains a PAZ domain, found in proteins involved in RNA-mediated gene silencing, and a C-terminal Piwi domain. The encoded protein is thought to function in maintenance of germline cells. Multiple transcript variants encoding different isoforms have been found for this gene. |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |              |
|-------|--------------|
| 基因ID  | 440822       |
| 基因名   | PIWIL3       |
| Swiss | Q7Z3Z3       |
| 别名    | PIWIL3;HIWI3 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)