

U2AF1 Rabbit pAb

货号: B19444

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200 IP: 1:500 - 1:1000
理论分子量	8kDa/19kDa/27kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HeLa,K-562,LO2,Mouse spleen,Mouse kidney,Rat spleen
细胞定位	Nucleus,Nucleus speckle
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human U2AF1 (NP_006749.1).
序列	MAEYLASIFGTEKDKVNCsfYFKIGACRHGDRCsRLHNKPTFSQTIALlNIYRNPQNSSQSADGLRCAVSDVEMQEHyDEF FEEVFTEMEEKYGEVEEMN

靶点信息

研究背景	This gene belongs to the splicing factor SR family of genes. U2 auxiliary factor, comprising a large and a small subunit, is a non-snRNP protein required for the binding of U2 snRNP to the pre-mRNA branch site. This gene encodes the small subunit which plays a critical role in both constitutive and enhancer-dependent RNA splicing by directly mediating interactions between the large subunit and proteins bound to the enhancers. Alternatively spliced transcript variants encoding different isoforms have been identified.
------	---

基因ID	7307
基因名	
Swiss	Q01081
别名	U2AF1;FP793;RN;RNU2AF1;U2AF35;U2AFBP

产品验证



Western blot analysis of U2AF1 expressed in Mouse spleen,Rat spleen,293T using U2AF1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn