

RAVER1 Rabbit pAb

货号: B19456

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	8kDa/63kDa/77kDa
实测分子量	64kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HeLa
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 390-490 of human RAVE R1 (NP_597709.2).
序列	GPALSTALLQLALQTQGQKKPGILGDSPLGALQPGAQPANPLLGELPAGGGLPPELPPIRRGKPPPLLPSVLGPAGGDREAL GLGPPAAQLTPPPAPVGLRG

靶点信息

研究背景	Cooperates with PTBP1 to modulate regulated alternative splicing events. Promotes exon skipping. Coope rates with PTBP1 to modulate switching between mutually exclusive exons during maturation of the TPM1 pre-mRNA (By similarity).
基因ID	125950

基因名	RAVER1
Swiss	Q8IY67
别名	RAVER1

产品验证



Western blot analysis of RAVR1 expressed in 293T,HeLa using RAVR1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn