

RCC2 Rabbit pAb

货号: B19467

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	56kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Jurkat
细胞定位	Chromosome,Cytoplasm,Nucleus,centromere,cytoskeleton,nucleolus,spindle
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 323-522 of human RCC2 (NP_001129676.1).
序列	QILPVPNVVVRDVACGANHTLVLDSQKRVFSWGGYGRLGHAEQKDEMPRLVKLDFPGRGASQIYAGYTCSFAVSEVGGLFFWGATNTSRESTMYPKAVQDLCGWIRISLACGKSSIIVAADESTISWGPSPTFGELGYGDHKPKSSTAAQEVKTLDGIFSEQVAMGYSHSLVIARDESETEKEKIKKLPEYNPRTL

靶点信息

研究背景	The protein encoded by this gene is a guanine exchange factor that is active on RalA, a small GTPase. The encoded protein and RalA are both essential for proper kinetochore-microtubule function in early mitosis. This protein has been shown to be a biomarker for colorectal cancer.
基因ID	55920

基因名	RCC2
Swiss	Q9P258
别名	RCC2;TD-60

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn