

NLN Rabbit pAb

货号: **B19471**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	80kDa
实测分子量	80kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,Rat liver
细胞定位	Cytoplasm,Mitochondrion intermembrane space
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 38-260 of human NLN (NP_065777.1).
序列	SSYTVAGRNVLRWDLSPEDIKTRTEELIVQTKQVYDAVGMLGIEEVTYENCLQALADVEVKYIVERTMLDFPQHVSSDKEV RAASTEADKRLSRFDIEMSMRGDIFERIVHLQETCDLGKIKPEARRYLESIKMGKRNGHLHLPQVQNEIKSMKKRMSEL DFNKNLNEDDTFLVFSKAELGALPDDFIDSLEKTDDDKYKITLKYPHYFPVMKKCCIPET

靶点信息

研究背景	This gene encodes a member of the metallopeptidase M3 protein family that cleaves neurotensin at the P ro10-Tyr11 bond, leading to the formation of neurotensin(1-10) and neurotensin(11-13). The encoded pro tein is likely involved in the termination of the neurotensinergic signal in the central nervous system and i n the gastrointestinal tract.
基因ID	57486

基因名	NLN
Swiss	Q9BYT8
别名	NLN;AGTBP;EP24.16;MEP;MOP

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn