

TRAFD1 Rabbit pAb

货号: B19487

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	19kDa/64kDa
实测分子量	78kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A375,HeLa
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human TRAFD1 (NP_006691.1).
序列	MAEFLDDQETRLCDNCKKEIPVFNFTHIHCQRNIGMCPTCKEFPKSDMETHMAAEHCQVTCKCNKKLEKRLKKHEETE CPLRLAVCQHCDLELSILKLKEHEDYCGARTELCGNCGRNVLVKDLKTHPEVCGREGEEKRNEVAIPPNAVYDESWGQDGI WIASQLLRQIEALDPPMRLPRRPLRAFESDVFHNRRTTNQRNITAQVSIQNNLFEEQERQERNRGQQPPKEGGEESEANLDF MLALSLQNEGQASSVAEQ

靶点信息

研究背景	The innate immune system confers host defense against viral and microbial infection, and TRAFD1 is a negative feedback regulator that controls excessive immune responses (Sanada et al., 2008 [PubMed 18849341]).[supplied by OMIM, Dec 2009]
基因ID	10906

基因名	TRAFD1
Swiss	O14545
别名	TRAFD1;FLN29

产品验证



Western blot analysis of TRAFD1 expressed in A375,HeLa using TRAFD1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn