

PITPNM1 Rabbit pAb

货号: B19491

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	134kDa
实测分子量	135kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,Mouse brain,Mouse eye,Rat brain
细胞定位	Cleavage furrow,Cytoplasm,Endoplasmic reticulum membrane,Golgi apparatus,Golgi stack membrane,Lipid droplet,Midbody,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 110-190 of human PITPNM1 (NP_004901.2).
序列	EIETYYLPDGGQQPNVFNLGAERRQRILDTIDIVRDAVAPGEYKAEDPRLYHSVKTRGRPLSDDWARTAAQTGPLMCAY

靶点信息

研究背景	PITPNM1 belongs to a family of membrane-associated phosphatidylinositol transfer domain-containing proteins that share homology with the Drosophila retinal degeneration B (rdgB) protein (Ocaka et al., 2005 [PubMed 15627748]).[supplied by OMIM, Mar 2008]
基因ID	9600

基因名	PITPNM1
Swiss	O00562
别名	PITPNM1;DRES9;NIR2;PITPNM;RDGB;RDGB1;RDGBA;RDGBA1;Rd9

产品验证



Western blot analysis of PITPNM1 expressed in Mouse brain,Rat brain,LO2 using PITPNM1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn