

IFT88 Rabbit pAb

货号: B19494

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	92kDa/93kDa/94kDa
实测分子量	108kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,BxPC-3,Mouse brain,Mouse eye
细胞定位	Cytoplasm,centriole,centrosome,cilium basal body,cytoskeleton,microtubule organizing center
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 534-833 of human IFT8 8 (NP_783195.2).
序列	GLTYEKLNRLEALDCFLKLHAILRNSAEVLYQIANIYELMENPSQAIEWLMQVVSVIPTDPQVLSKLGELYDREGDKSQAF QYYYESYRYFPCNIEVIEWLGAYYIDTQFWEKAIQYFERASLIQPTQVKWQLMVASCFRRSGNYQKALDITYKDTHRKFPEN VECLRFLVRLCTDLGLKDAQEYARKLKRLEKMKEIREQRIKSGRDGSGGSRGKREGSASGDSGQNYSSASSKGERLSARLR ALPGTNEPYESSNKEIDASYVDPLGPQIERPKTAAKKRIDEEDFADEELGDDLLPE

靶点信息

研究背景	This gene encodes a member of the tetratrico peptide repeat (TPR) family. Mutations of a similar gene in mouse can cause polycystic kidney disease. Two transcript variants encoding distinct isoforms have been identified for this gene.
------	--

基因ID	8100
基因名	IFT88
Swiss	Q13099
别名	IFT88;D13S1056E;DAF19;TG737;TTC10;hTg737

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn