

WSB2 Rabbit pAb

货号: B19520

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	21kDa/45kDa/47kDa
实测分子量	40-43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Rat kidney,Rat spleen
细胞定位	cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human WSB2 (NP_061109.1).
序列	MEAGEEPLLLAELKGRPHQFDWKSSCETWSVAFSPDGSWFAWSQGHCVKLIPWPLEEQFIPKGFEAKSRSSKNETKG RGSPKEKTLDCGQIVWGLAFSPWPSPSRKLWARHHPQVPDVSVCLVLATGLNDGQIKIWEVQTGLLLNLSGHQDVVRD LSFTPSGSLILVSASRDKTLRIWDLNKHGKQIQVLSGHLQWV

靶点信息

研究背景	This gene encodes a member of the WD-protein subfamily. The encoded protein contains five WD-repeats spanning most of the protein and an SOCS box in the C-terminus. The SOCS box may act as a bridge between specific substrate-binding domains and E3 ubiquitin protein ligases. Three transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	55884
基因名	WSB2
Swiss	Q9NYS7
别名	WSB2;SBA2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn