

ST6GALNAC4 Rabbit pAb

货号: B19560

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:2000
理论分子量	34kDa
实测分子量	34kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Mouse liver
细胞定位	Golgi apparatus membrane,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 120-210 of human ST6GALNAC4 (NP_778204.1).
序列	TLRVVSHTSVPLLLRNYSHYFQKARDTLYMWVGQGRHMDRVLGGRTYRTLQLTRMYPGLQVYTFTERMMAYCDQIFQDETGKNRRQSGSF

靶点信息

研究背景	The protein encoded by this gene is a type II membrane protein that catalyzes the transfer of sialic acid from CMP-sialic acid to galactose-containing substrates. The encoded protein prefers glycoproteins rather than glycolipids as substrates and shows restricted substrate specificity, utilizing only the trisaccharide sequence Neu5Ac-alpha-2,3-Gal-beta-1,3-GalNAc. In addition, it is involved in the synthesis of ganglioside G _{D1A} from GM _{1B} . The encoded protein is normally found in the Golgi apparatus but can be proteolytically processed to a soluble form. This protein is a member of glycosyltransferase family 29. Transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	27090
基因名	ST6GALNAC4
Swiss	Q9H4F1
别名	ST6GALNAC4;IV;SIAT3-C;SIAT3C;SIAT7-D;SIAT7D;ST6GALNACIV;ST6GalNAc

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn