

TP53I11 Rabbit pAb

货号: B19568

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	21kDa
实测分子量	21KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse heart
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human TP53I11 (NP_001245249.1).
序列	MAAKQPPPLMKKHSQTDLVSRKTRKILGVGGEDDDGEVHRSKISQVLGNEIKFTIREPLGLRVWQFVSA

靶点信息

研究背景	Predicted to be involved in negative regulation of cell population proliferation. Predicted to be integral component of membrane.
基因ID	9537

基因名	TP53I11
Swiss	O14683
别名	TP53I11;PIG11

产品验证



Western blot analysis of TP53I11 expressed in Mouse brain,Mouse heart using TP53I11 Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Bloc king buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn