

EPB41L3 Rabbit pAb

货号: B19614

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	69kDa/96kDa/98kDa/120kDa
实测分子量	121kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Mouse kidney
细胞定位	Cell junction,Cell membrane,Cytoplasm,Cytoplasmic side,Peripheral membrane protein,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 656-865 of human EPB4 1L3 (NP_001268463.1).
序列	GDAAAQPAFTGIKGKEGSALTEGAKEEGGEEVAKAVLEQEETAAASRERQEEQSAAIHISETLEQKPHFESSTVKTETISFG SVSPGGVKLEISTKEVPVVHTETKTITYESSQVDPGTDLEPGVLSAQTTSETTSTTTTTHITKTVKGGISETRIEKRIVITGD ADIDHDQALAQAIKEAKEQHPDMSVTKVVVHKETEITPEDGED

靶点信息

研究背景	Tumor suppressor that inhibits cell proliferation and promotes apoptosis. Modulates the activity of protein arginine N-methyltransferases, including PRMT3 and PRMT5.
基因ID	23136

基因名	EPB41L3
Swiss	Q9Y2J2
别名	EPB41L3;4.1B;DAL-1;DAL1

产品验证



Western blot analysis of EPB41L3 expressed in U-87MG, Mouse kidney using EPB41L3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn