

TRIM45 Rabbit pAb

货号: B19615

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IF/ICC
推荐浓度	IF/ICC: 1:50 - 1:200
理论分子量	62kDa/64kDa
实测分子量	64KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 100-400 of human TRIM 45 (NP_001139107.1).
序列	QVDLPMGGVKALTIDHLAVNDVMLES LRGEGQGLVCDLCNDREVEKRCQTCKANLCHFCCQAHRQKKTYYHTMVDLKDLKGYSRIGKPILCPVHPAEELRLFCEFCDRPVCQDCVVGEHREHPCDFTSNVIHKHGDSVWELLKGTQPHVEALEEALAQIHIINSALQKRVEAVAADVRTFSEGYIKAIEEHRDKLLKQLEDIRAQKENS LQLQKAQLEQLLADMRTGVEFTEHLLTSGSDLEILITKRVVVERLRKLNKVQYSTRPGVNDKIRFCPQEKAGQCRGYEIDLHRAREKQTA

靶点信息

研究背景	This gene encodes a member of the tripartite motif family. The encoded protein may function as a transcriptional repressor of the mitogen-activated protein kinase pathway. Alternatively spliced transcript variants have been described.
基因ID	80263

基因名	TRIM45
Swiss	Q9H8W5
别名	TRIM45;RNF99

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn