

# NUCB2 Rabbit pAb

货号: B19639

产品信息

|       |                                                                                                                |
|-------|----------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Rat                                                                                                      |
| 宿主    | Rabbit                                                                                                         |
| 克隆性   | Polyclonal                                                                                                     |
| 预测反应  |                                                                                                                |
| 应用    | WB IHC                                                                                                         |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:100 - 1:200                                                         |
| 理论分子量 | 50kDa                                                                                                          |
| 实测分子量 | 50kDa                                                                                                          |
| 形式    | Liquid                                                                                                         |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.           |
| 偶联物   | Unconjugated                                                                                                   |
| 阳性对照  | LO2,HeLa,Jurkat,SGC-7901                                                                                       |
| 细胞定位  | Cytoplasm,Endoplasmic reticulum,Golgi apparatus,Membrane,Nucleus envelope,Peripheral membrane protein,Secreted |
| 纯化    | Affinity purification                                                                                          |

抗原信息


|      |                                                                                                                     |
|------|---------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 321-420 of human NUCB2 (NP_005004.1). |
| 序列   | ATEKKEFLEPDSWETLDQQQFFTEELKEYENIIALQENELKKKADELQKQKEELQRQHDQLEAQKLEYHQVIQQMEQKKLQQGIPPSGPAGELKFEPHI                 |

靶点信息

|      |                                                                                                                                                                                                                                                     |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a protein with a suggested role in calcium level maintenance, eating regulation in the hypothalamus, and release of tumor necrosis factor from vascular endothelial cells. This protein binds calcium and has EF-folding domains. |
| 基因ID | 4925                                                                                                                                                                                                                                                |

|       |                      |
|-------|----------------------|
| 基因名   | NUCB2                |
| Swiss | P80303               |
| 别名    | NUCB2;HEL-S-109;NEFA |

产品验证



Western blot analysis of NUCB2 expressed in LO2,HeLa,Jurkat using NUCB2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)