

RASSF5 Rabbit pAb

货号: B19651

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	30kDa/37kDa/47kDa
实测分子量	36-52kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	THP-1,A-549,293T,Jurkat,Mouse brain
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-265 of human RASSF5 (NP_872606.1).
序列	MTVDSSMSSGYCSLDEELEDCEFTAKTTFFRNAQSKHLSKNVCKPVEETQRPPTLQEIKQKIDSYNTRKNCLGMKLS EDGTYTGFIKVHLKLRRPVTVPAGIRPQSIYDAIKEVNLAATTDKRTSFYLPDAIKQLHISSTTTVSEVIQGLLKKFMVVDNPQK FALFKRIHKDGGVLFQKLSIADRPLYRLLAGPDTEVLSFVLKENETGEVEWDAFSIPELQNF LTILEKEEQDKIQVQKKYD KFRQKLEEALRESQGKPG

靶点信息

研究背景	This gene is a member of the Ras association domain family. It functions as a tumor suppressor, and is in activated in a variety of cancers. The encoded protein localizes to centrosomes and microtubules, and as sociates with the GTP-activated forms of Ras, Rap1, and several other Ras-like small GTPases. The protei n regulates lymphocyte adhesion and suppresses cell growth in response to activated Rap1 or Ras. Multip le transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	83593
基因名	RASSF5
Swiss	Q8WWW0
别名	RASSF5;Maxp1;NORE1;NORE1A;NORE1B;RAPL;RASSF3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn