

BIRC7 Rabbit pAb

货号: B19652

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:100
理论分子量	30kDa/32kDa
实测分子量	35kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Mouse skeletal muscle
细胞定位	Cytoplasm,Golgi apparatus,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-265 of human BIRC7 (NP_647478.1).
序列	MGPKDSAKCLHRGPQPSHWAAGDGPTQERCGRPRLGSPVLGLDTCRAWDHVDGQILGQLRPLTEEEEEEGAGATLSRG PAFPGMGSEELRLASFYDWPLTAEVPPELLAAAGFFHTGHQDKVRCFFCYGGLQSWKRGDDPWTEHAKWFPSCQFLLR SKGRDFVHSVQETHSQLLGSWDPWEEPDAAPVAPSVPASGYPELPTPRREVQSESAQEPGGVSPAQAQRAWWVLEPP GARDVEAQLRRLQEERTCKVCLDRAVSIVFV

靶点信息

研究背景	This gene encodes a member of the inhibitor of apoptosis protein (IAP) family, and contains a single copy of a baculovirus IAP repeat (BIR) as well as a RING-type zinc finger domain. The BIR domain is essential for inhibitory activity and interacts with caspases, while the RING finger domain sometimes enhances anti apoptotic activity but does not inhibit apoptosis alone. Elevated levels of the encoded protein may be associated with cancer progression and play a role in chemotherapy sensitivity. Alternative splicing results in multiple transcript variants
基因ID	79444
基因名	BIRC7
Swiss	Q96CA5
别名	BIRC7;KIAP;LIVIN;ML-IAP;MLIAP;RNF50

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn