

HIRA Rabbit pAb

货号: B19681

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	89kDa/111kDa
实测分子量	112kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,HeLa,Mouse kidney
细胞定位	Nucleus,PML body
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 250-550 of human HIRA (NP_003316.3).
序列	AQIIEREGWKTNMDFVGHRKAVTVVKFNPKIFKKKQKNGSSAKPSCPYCCCAVGSKDRSLSVWLTCLKRPLVVIHELFDKS IMDISWTLNGLGILVCSMDGSVAFLDFSQDELGDPLSEEEKSRIHQSTYGKSLAIMTEAQLSTAVIENPEMLKYQRRQQQQ QLDQKSAATREMGSATSVAGVVNGESLEDIRKNLLKKQVETRTADGRRRRITPLCIAQLDTGDFSTAFFNSIPLSGSLAGTM LSSHSPQLPLDSSTPNSFGASKPCTEPVVAASARPAGDSVNKDSMNATSTPAALSP

靶点信息

研究背景	This gene encodes a histone chaperone that preferentially places the variant histone H3.3 in nucleosome s. Orthologs of this gene in yeast, flies, and plants are necessary for the formation of transcriptionally silent heterochromatin. This gene plays an important role in the formation of the senescence-associated heterochromatin foci. These foci likely mediate the irreversible cell cycle changes that occur in senescent cells . It is considered the primary candidate gene in some haploinsufficiency syndromes such as DiGeorge syndrome, and insufficient production of the gene may disrupt normal embryonic development.
------	--

基因ID	7290
基因名	HIRA
Swiss	P54198
别名	HIRA;DGCR1;TUP1;TUPLE1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn