

RBP4 Rabbit pAb

货号: B19689

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	23kDa
实测分子量	19kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620,Mouse liver
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 19-201 of human RBP4 (NP_006735.2).
序列	ERDCRVSSFRVKENFDKARFSGTWYAMAKKDPEGLFLQDNIVAEFVSDETGQMSATAKGRVRLNNWVDCADMVGTF TDTEPAKFKMKYWGVASFLLQKGNDDHWIVDTDYDTYAVQYSCRLLNLDGTCADSYSFVFSRDPNGLPPEAQKIVRQR QEELCLARQYRLIVHNGYCDGRSERNLL

靶点信息

研究背景	This protein belongs to the lipocalin family and is the specific carrier for retinol (vitamin A alcohol) in the blood. It delivers retinol from the liver stores to the peripheral tissues. In plasma, the RBP-retinol complex interacts with transthyretin which prevents its loss by filtration through the kidney glomeruli. A deficiency of vitamin A blocks secretion of the binding protein posttranslationally and results in defective delivery and supply to the epidermal cells.
------	--

基因ID	5950
基因名	RBP4
Swiss	P02753
别名	RBP4;MCOPCB10;RDCCAS

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn