

Alpha-1 Antitrypsin (SERPINA1) Rabbit pAb

货号: B19695

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	34kDa/40kDa/46kDa
实测分子量	45kDa,50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,HT-29,Mouse lung
细胞定位	Endoplasmic reticulum,Secreted,extracellular matrix,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 25-315 of human Alpha-1 Antitrypsin (Alpha-1 Antitrypsin (SERPINA1)) (NP_000286.3).
序列	EDPQGDAAQKTDTSHTDQDHPTFNKITPNLAFAFSLYRQLAHQSNSTNIFFSPVSIATAFAMLSLGTKADTHDEILEGLN FNLTEIPEAQIHEGFQELLRTLNPDSQLQLTTGNGLFLSEGLKLVDFLEDVKKLYHSEAFTVNFGDTEEAKKQINDYVEK GTQGKIVDLVKELDRDTVFALVNYIFFKGKWERPFVKTDEEEDFHVDQVTTVKVPMMKRLGMFNIQHCKKLSSWVLLM KYLGNATAIFFLPDEGKLQHLENELTHDIITKFLNEDRRSASLHLPKL

靶点信息

研究背景	The protein encoded by this gene is secreted and is a serine protease inhibitor whose targets include elas tase, plasmin, thrombin, trypsin, chymotrypsin, and plasminogen activator. Defects in this gene can caus e emphysema or liver disease. Several transcript variants encoding the same protein have been found for this gene.
------	--

基因ID	5265
基因名	SERPINA1
Swiss	P01009
别名	SERPINA1;A1A;A1AT;AAT;PI;PI1;PRO2275;alpha1AT

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn