

PALLD Rabbit pAb

货号: B19743

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IP IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IP: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	42kDa/58kDa/73kDa/85-150kDa
实测分子量	150kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549
细胞定位	Cell junction,Cell projection,Cytoplasm,Z line,cytoskeleton,focal adhesion,lamellipodium,myofibril,ruffle,s arcomere
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 383-672 of human PALL D (NP_001159582.1).
序列	AVNQRGRSPRSPSGHPHVRRPRSRSDSGDENIQERFFRPHFLQAPGDLTVQEGKLCRMDCKVSGLPTPDLWSWQLD GKPVRPDSAHKMLVRENGVHSLIIEPVTSRDAGIYTCIATNRAGQNSFSLELVVAAKEAHKPPVFIEKLQNTGVADGYPVRL ECRVLGVPQQIFWKKENESLTHSTDRVSMHQDNHGYICLLIQGATKEDAGWYTVSAKNEAGIVSCTARLDVYTQWHQQ SQSTKPKKVRPSASRYAALSDQGLDIKAAFQPEANPSHLTLNTALVEEDL

靶点信息

研究背景	This gene encodes a cytoskeletal protein that is required for organizing the actin cytoskeleton. The protein is a component of actin-containing microfilaments, and it is involved in the control of cell shape, adhesion, and contraction. Polymorphisms in this gene are associated with a susceptibility to pancreatic cancer type 1, and also with a risk for myocardial infarction. Alternative splicing results in multiple transcript variants.
基因ID	23022
基因名	PALLD
Swiss	Q8WX93
别名	PALLD;CGI-151;CGI151;MYN;PNCA1;SIH002;palladin

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn