

ARAP2 Rabbit pAb

货号: B19769

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	193kDa
实测分子量	170kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,K-562,Mouse brain
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1405-1704 of human ARAP2 (NP_056045.2).
序列	LAEPGSAYLVVKRFLTADTIKHCSRSTLGSIKEGILKIKEEPSKILSGNKFQDRYFVLRDGFLLYKDVKSSKHKDKMFSLS MKFYRGVKKMKPPTSWGLTAYSEKHHWHLCCDSSRTQTEWMTSIFIAQHEYDIWPPAGKERKRSITKNPKIGGLPLIQ HEGNATLARKNIESARAELERLRLSEKCDKESVDSSLKERASMAHLCLEHKDDKLRNRPRKHSFNCLEDTEPEAPLGQPK GHKGLKTLRKTEDRNSKATLDSHKLPSRVIEELNVVLQRSRTLPELQDEQILK

靶点信息

研究背景	The protein encoded by this gene contains ARF-GAP, RHO-GAP, ankyrin repeat, RAS-associating, and pleckstrin homology domains. The protein is a phosphatidylinositol (3,4,5)-trisphosphate-dependent Arf6 GAP that binds RhoA-GTP, but it lacks the predicted catalytic arginine in the RHO-GAP domain and does not have RHO-GAP activity. The protein associates with focal adhesions and functions downstream of RhoA to regulate focal adhesion dynamics.
------	---

基因ID	116984
基因名	ARAP2
Swiss	Q8WZ64
别名	ARAP2;CENTD1;PARX

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn