

UBR4 Rabbit pAb

货号: B19797

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	23kDa/270kDa/571kDa/573kDa/575kDa
实测分子量	600kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Cytoplasm,Membrane,Multi-pass membrane protein,Nucleus,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 4884-5183 of human UBR4 (NP_065816.2).
序列	NIVHYDCHLAAVRLARGREEWESAALQNANTKCNGLLPVWGPHVPESAFATCLARHNTYLQECTGQREPTYQLNIHDIKLLFLRFAMEQSFSADTGGGGRESNIHLIPYIIHTVLYVLNLTTRATSREEKNLQGFLEQPKEKWVESAFEVDGPYYFTVLALHILPPEQWRATRVEILRRLLVTSQARAVAPGGATRLTDKAVKDYSAYRSSLLFWALVDLIYNMFKKVPTSNTEGGWSCSLAEYIRHNDMPIYEAADKALKTFQEEFMPVETFSEFLDVAGLLSEITDPESFLKDLLNSVP

靶点信息

研究背景	The protein encoded by this gene is an E3 ubiquitin-protein ligase that interacts with the retinoblastoma-associated protein in the nucleus and with calcium-bound calmodulin in the cytoplasm. The encoded protein appears to be a cytoskeletal component in the cytoplasm and part of the chromatin scaffold in the nucleus. In addition, this protein is a target of the human papillomavirus type 16 E7 oncoprotein.
------	--

基因ID	23352
基因名	UBR4
Swiss	Q5T4S7
别名	UBR4;RBAF600;ZUBR1;p600

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn