

PIGT Rabbit pAb

货号: B19824

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	41kDa/43kDa/54kDa/57kDa/59kDa/65kDa
实测分子量	66kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A-431,LO2,B cells,Mouse brain,Mouse testis,Rat brain
细胞定位	Endoplasmic reticulum membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 260-460 of human PIGT (NP_057021.2).
序列	LFRMFSRTLTEPCPLASESRVYVDITTYNQDNETLEVHPPPTTTYQDVILGTRKTYAIYDLLDTAMINNSRNLNIQLKWKRPPE NEAPPVPFLHAQRYVSGYGLQKGELSTLLYNTHPYRAFPVLLDTPWYLRLYVHTLTITSKGKENKPSYIHYQPAQDRLQP HLEMLIQLPANSVTKVSIQFERALLKWTEYTPDP

靶点信息

研究背景	This gene encodes a protein that is involved in glycosylphosphatidylinositol (GPI)-anchor biosynthesis. The GPI-anchor is a glycolipid found on many blood cells and serves to anchor proteins to the cell surface. This protein is an essential component of the multisubunit enzyme, GPI transamidase. GPI transamidase mediates GPI anchoring in the endoplasmic reticulum, by catalyzing the transfer of fully assembled GPI units to proteins. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
------	--

基因ID	51604
基因名	PIGT
Swiss	Q969N2
别名	PIGT;CGI-06;MCAHS3;NDAP;PNH2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn