

PAI-1/Serpin E1 Rabbit pAb

货号: **B19829**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	43kDa/45kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-1080,SH-SY5Y,Mouse liver,Rat liver
细胞定位	Secreted
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 24-158 of human PAI-1/ Serpin E1 (NP_000593.1).
序列	VHHPPSYVAHLASDFGVRVFQQVAQASKDRNVVFSYPGVASVLAMLQLTTGGETQQQIQAAAMGFKIDDKGMAPALRHL YKELMGPWNKDEISTTDAIFVQRDLKLVQGFMPIFFRLFRSTVKQVDFSEVERARFI

靶点信息

研究背景	This gene encodes a member of the serine proteinase inhibitor (serpin) superfamily. This member is the p rincipal inhibitor of tissue plasminogen activator (tPA) and urokinase (uPA), and hence is an inhibitor of fib rinolysis. Defects in this gene are the cause of plasminogen activator inhibitor-1 deficiency (PAI-1 deficien cy), and high concentrations of the gene product are associated with thrombophilia. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	5054
基因名	SERPINE1
Swiss	P05121
别名	PAI;PAI-1;PAI1;PLANH1;SERPINE1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn