

PKLR Rabbit pAb

货号: B19834

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	58kDa/61kDa
实测分子量	62kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,HepG2,293T,Mouse kidney,Mouse liver,Mouse heart,Rat kidney
细胞定位	cytoplasm,cytosol,extracellular exosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 446-556 of human PKLR (NP_000289.1).
序列	PLSRDPTEVTAIGAVEAAFKCCAAAIIVLTTTGRSAQLLSRYRPRAAVIAVTRSAQAARQVHLCRGVFPLLYREPPEAIWADDVDRRVQFGIESGKLRGFLRVGDLVIVVT

靶点信息

研究背景	The protein encoded by this gene is a pyruvate kinase that catalyzes the transphosphorylation of phohsp hoenolpyruvate into pyruvate and ATP, which is the rate-limiting step of glycolysis. Defects in this enzym e, due to gene mutations or genetic variations, are the common cause of chronic hereditary nonspherocy tic hemolytic anemia (CNSHA or HNSHA). Multiple transcript variants encoding different isoforms have be en found for this gene.
基因ID	5313

基因名	PKLR
Swiss	P30613
别名	PKLR;PK1;PKL;PKR;PKRL;RPK

产品验证



Western blot analysis of PKLR expressed in Mouse kidney,Rat kidney,BxPC-3 using PKLR Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn