

Septin 11 Rabbit pAb

货号: B19848

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:2000 - 1:6000
理论分子量	49kDa/50kDa
实测分子量	49KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HeLa,Jurkat,NIH/3T3,Mouse brain,Mouse kidney,Rat brain,Rat lung,Rat spleen
细胞定位	Cell junction,Cell projection,Cytoplasm,axon,cytoskeleton,dendritic spine,synapse
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human Septin 11 (NP_060713.1).
序列	MAVAVGRPSNEELRNLSLSGHVGFDSLDPQLV NKSTSQGFCFNILCVGETGIGKSTLMDTLFNTKFESDPATHNEPGVRLKARSYELQESNVRLKLTIVDTVGF GDQINKDDSYKPIVEY

靶点信息

研究背景	SEPT11 belongs to the conserved septin family of filament-forming cytoskeletal GTPases that are involved in a variety of cellular functions including cytokinesis and vesicle trafficking (Hanai et al., 2004 [PubMed 15196925]; Nagata et al., 2004 [PubMed 15485874]).[supplied by OMIM, Jul 2009]
基因ID	55752

基因名	SEPT11
Swiss	Q9NVA2
别名	SEPT11;septin-11

产品验证



Western blot analysis of Septin 11 expressed in Mouse brain,Rat brain,293T using Septin 11 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120 S.

实验步骤

访问官网浏览详情: www.ablybio.cn