

VISTA Rabbit pAb

货号: B19858

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	33kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,LO2,A-549
细胞定位	Cell membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 33-194 of human VISTA (NP_071436.1).
序列	FKVATPYSLYVCPegQNVTLTCRLLGPVDKGHDVTFYKTWYRSSRGEVQTCSERRPIRNLTFQDLHLHHGGHQAANTSH DLAQRHGLESASDHHGNFSITMRNLTLDSGLYCCLVVEIRHHHSEHRVHGAMELQVQTGKDAPSNVCVYPSSSQDSEN ITAA

靶点信息

研究背景	V-type immunoglobulin domain-containing suppressor of T-cell activation(VISTA) also known as platelet re ceptor Gi24 , stress-induced secreted protein-1 (Sisp-1) and B7
基因ID	64115

基因名	VSIR
Swiss	Q9H7M9
别名	VSIR;B7-H5;B7H5;C10orf54;DD1alpha;GI24;PD-1H;PP2135;SISP1;VISTA

产品验证



Western blot analysis of VISTA expressed in HeLa,LO2,A-549 using VISTA Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn