

STX2 Rabbit pAb

货号: B19871

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	32kDa/33kDa
实测分子量	36kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,DU145,SKOV3,A-549,Mouse heart,Mouse testis,Rat liver,Rat heart,Rat testis
细胞定位	Membrane,Single-pass type IV membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human STX2 (NP_001971.2).
序列	MRDRLPDLTACRKNDDGDTVVVVEKDHFMDDFFHQVEEIRNSIDKITQYVEEVKKNHSIILSAPNPEGKIKEELEDLNKEIK KTANKIRAKLKAIEQSFDQDESGNRTSVDLRIRRTQHSVLSRKFEAMA EYNEAQT LFRERSKGRIQRQLEITGR TTTDDEL EEMLESGKPSIFTSDIISDSQITRQALNEIESRHKDIMKLETSIRELHEMFMDMAMFVETQGEMINNIERNVMNATDYVEHAK EETKKAIKYQSKARRKL

靶点信息

研究背景	The product of this gene belongs to the syntaxin/epimorphin family of proteins. The syntaxins are a large protein family implicated in the targeting and fusion of intracellular transport vesicles. The product of this gene regulates epithelial-mesenchymal interactions and epithelial cell morphogenesis and activation. Alternatively spliced transcript variants encoding different isoforms have been identified.
------	---

基因ID	2054
基因名	STX2
Swiss	P32856
别名	STX2;EPIM;EPM;STX2A;STX2B;STX2C

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn