

PCIF1 Rabbit pAb

货号: B19888

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:1000 - 1:5000 IHC: 1:50 - 1:200 IP: 1:50 - 1:200
理论分子量	81kDa
实测分子量	80KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	intercellular bridge,microtubule cytoskeleton,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PCIF1 (NP_071387.1).
序列	MANENHGSPREEASLLSHSPGTSNQSQPCSPKPIRLVQDLPEELVHAGWEKCWSRRENRPYYFNRFTNQSLWEMPVLGQHDVISDPLGLNATPLPQDSSLVETPPAENKPRKRQLSEEQPSGNGVKKPKIEIPVPTGTGQSVSPSSIPGTPTLKMWGTSPEDKQQAALLRPTEVYWDLDIQTNAVIKHRGPSEVLPPHPE

靶点信息

研究背景	Cap-specific adenosine methyltransferase that catalyzes formation of N(6,2'-O-dimethyladenosine cap (m6A)m by methylating the adenosine at the second transcribed position of capped mRNAs. Recruited to the early elongation complex of RNA polymerase II (RNAPII via interaction with POLR2A and mediates formation of m6A)m co-transcriptionally.
------	--

基因ID	63935
基因名	PCIF1
Swiss	Q9H4Z3
别名	CAPAM;hCAPAM;hPCIF1;C20orf67;PPP1R121;PCIF1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn