

MRPS9 Rabbit pAb

货号: B19907

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
理论分子量	45kDa
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,U-87MG,A-431,HT-29
细胞定位	Mitochondrion
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 237-396 of human MRP S9 (NP_872578.1).
序列	EEEFVQRFRRSVTLSEKKQLIEPVQYDEQGMAFSKSEGKRKTAKAEIVYKHGSGRIKVNIDYQLYFPITQDREQLMFPFH FVDRLGKHDVTCTVSGGGRSAQAGAIRLAMAKALCSFVTEDEVWMRQAGLLTTPRVRRERKKPGQEGARRKFTWKKR

靶点信息

研究背景	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein.
基因ID	64965
基因名	MRPS9
Swiss	P82933
别名	MRPS9;MRP-S9;RPMS9;S9mt

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn