

CHAMP1 Rabbit pAb

货号: B19918

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	89kDa
实测分子量	130kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A-431,Jurkat
细胞定位	Chromosome,Cytoplasm,Nucleus,centromere,cytoskeleton,kinetochore,spindle
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 563-812 of human CHAMP1 (NP_115812.1).
序列	LPKSALFSESQKAVELGDELQIDAIDDQKCDILVQEELLASPKKLEDTLFPSSKKLKKDNQESSDAELSSSEYIKTDLDAM DIKGQESSSDQEQQVDVESIDFSKENKMDMTSPEQSRNVLQFTEEKEAFISEEEIAKYMKRGKGKYYCKICCCRAMKKGAVL HHLVNHVHSPYKCTICGKAFLLKSLKNHVAAHGQSLKLCPRCNFESNFPRGFKKHLTHCQSRHNEEANKKLMEALEPP LEEQQI

靶点信息

研究背景	This gene encodes a zinc finger protein that functions as a regulator of chromosome segregation in mitosis. The encoded protein is required for correct alignment of chromosomes on the metaphase plate, and plays a role in maintaining the attachment of sister kinetochores to microtubules from opposite spindle poles. Mutations in this gene are associated with autosomal dominant mental retardation-40. [provided by RefSeq, Oct 2016]
基因ID	283489
基因名	CHAMP1
Swiss	Q96JM3
别名	CHAMP1;C13orf8;CAMP;CHAMP;MRD40;ZNF828

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn