

MC3R Rabbit pAb

货号: B19947

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	36kDa
实测分子量	36kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse liver,Mouse kidney,Rat brain
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human MC3R (NP_063941.3).
序列	MNASCCLPVQPTLPNGSEHLQAPFFSNQSSSAFCEQVFIKPEVFLSLGIVSLLLENILVILAVVRNGNLHSPMYFFLCSLAVA DMLVSVSNALETIMIAIVHSDYLTTFEDQFIQHMDNIF

靶点信息

研究背景	This gene encodes a G-protein-coupled receptor for melanocyte-stimulating hormone and adrenocorticotrophic hormone that is expressed in tissues other than the adrenal cortex and melanocytes. This gene maps to the same region as the locus for benign neonatal epilepsy. Mice deficient for this gene have increased fat mass despite decreased food intake, suggesting a role for this gene product in the regulation of energy homeostasis. Mutations in this gene are associated with a susceptibility to obesity in humans.
基因ID	4159

基因名	MC3R
Swiss	P41968
别名	MC3R;BMIQ9;MC3;MC3-R;OB20;OQTL

产品验证



Western blot analysis of MC3R expressed in Mouse brain,Rat brain,Mouse liver using MC3R Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn