

RSU1 Rabbit pAb

货号: B19949

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	25kDa/31kDa
实测分子量	31kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,U-87MG,B cells,SKOV3,Mouse heart,Mouse lung,Mouse kidney,Rat heart,Rat lung
细胞定位	cytosol,extracellular exosome,focal adhesion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-277 of human RSU1 (NP_036557.1).
序列	MSKSLKKLVEESREKNQPEVDMSDRGISNMLDVNGLFTLSHITQLVLSHNKLTMVPPNIAELKNLEVLNFFNNQIEELPTQIS SLQKLKHLNLMNRNLNLTLPGRFGSLPALEVLDTYNNLSENSLPGNFFYLTTLRALYLSDNDFEILPPDIGKLTQLQILSLRDN DLISLPKEIGELTQLKELHIQGNRLTVLPPELGNLDLTGQKQVFKAENNPWVTPIADQFQLGVSHVFYIRSETYKYLYGRHM QANPEPPKKNNDKSKKISRKPLAAKNR

靶点信息

研究背景	This gene encodes a protein that is involved in the Ras signal transduction pathway, growth inhibition, and nerve-growth factor induced differentiation processes, as determined in mouse and human cell line studies. In mouse, the encoded protein was initially isolated based on its ability to inhibit v-Ras transformation. Multiple alternatively spliced transcript variants for this gene have been reported; one of these variants was found only in glioma tumors.
------	---

基因ID	6251
基因名	RSU1
Swiss	Q15404
别名	RSU1;RSP-1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn