

ZNF444 Rabbit pAb

货号: B20013

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse lung,Mouse spleen
细胞定位	nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-327 of human ZNF444 (NP_060807.2).
序列	MEVAVPVKQEAEGALDSPWHRFRFRHLGDAPGPREALGLLRALCRDWLRPEVHTKEQMLELLVLEQFLSALPADTQAW VCSRQPQSGEEAVALLEELWGPAASPDGSSATRVPDVTQGPATGGKEDSGMIPLAGTAPGAEGPAPGDSQAVRPYK QEPSSPLAPGLPAFLAAPGTTSCPECCKTSLKPAHLLRHRQSHSGEKPHACPECCKAFRRKEHLRRHRDTHPGSPGSPG PALRPLPAREKPHACCECGKTFYWREHLVRHRKTHSGARPFACWECGKGFGRRREHVLRHQRIHGAAAASAQGAVAPGPD GGGPFPPWPLG

靶点信息

研究背景	This gene encodes a zinc finger protein which activates transcription of a scavenger receptor gene involved in the degradation of acetylated low density lipoprotein (Ac-LDL) (PMID: 11978792). This gene is located in a cluster of zinc finger genes on chromosome 19 at q13.4. A pseudogene of this gene is located on chromosome 15. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011]
基因ID	55311
基因名	ZNF444
Swiss	Q8N0Y2
别名	EZF-2;EZF2;ZSCAN17;ZNF444

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn