

TRMT61A Rabbit pAb

货号: B20050

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	31kDa
实测分子量	33KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PC-3,U-251MG,A-549,Mouse spleen,Rat testis
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-289 of human TRMT6 1A (NP_689520.2).
序列	MSFVAYEELIKEGDTAILSLGHGAMVAVRVQRGAQTQTRHGVLRHSDLIGRPFGSKVTCGRGGWVYVLHPTPELWTLNL PHRTQILYSTDIALITMMLELRPGSVVCESGTGSGSVSHAIIRTIAPTGHLTVEFHQQRAEKAREEFQEHRVGRWVTVRTQ DVCRSFGFVSHVADAVFLDIPSPWEAVGHAWDALKVEGGRFCFSFPCIEQVQRTCQALAARGFSELSTLEVLPQVYNVR TVSLPPPDLTGTGDPAGSDTSPFRSGTPMKEAVGHTGYLTFATKTPG

靶点信息

研究背景	Catalytic subunit of tRNA (adenine-N(1--methyltransferase, which catalyzes the formation of N(1-methyladenine at position 58 (m1A58 in initiator methionyl-tRNA. Catalytic subunit of mRNA N(1-methyltransferase complex, which mediates methylation of adenosine residues at the N(1 position of a small subset of mRNAs: N(1 methylation takes place in tRNA T-loop-like structures of mRNAs and is only present at low stoichiometries.
------	--

基因ID	115708
基因名	TRMT61A
Swiss	Q96FX7
别名	TRMT61A;C14orf172;GCD14;Gcd14p;TRM61;hTRM61

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn