

SEC14L2 Rabbit pAb

货号: B20053

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:1000 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	36kDa/44kDa/46kDa
实测分子量	46kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse liver,Mouse lung,Rat liver
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 60-270 of human SEC14 L2 (NP_001191133.1).
序列	GRKVETITIIYDCEGLGLKHLWKPAVEAYGEFLCMFEENYPETLKRLFVVKAPKLFVPVAYNLIKPFALSEDTRKKIMVLGANWKE VLLKHISPDQVPVEYGGTMTDPDGNPKCKSKINYGGDIPRKYYVRDQVKQQYEHSVQISRGSSSHQVEYEILFPGCVLRWQ FMSDGADVGFGLFKTKMGERQRAGEMTEVLPNQRYNSHLVPEDGTL

靶点信息

研究背景	This gene encodes a cytosolic protein which belongs to a family of lipid-binding proteins including Sec14p , alpha-tocopherol transfer protein, and cellular retinol-binding protein. The encoded protein stimulates squalene monooxygenase which is a downstream enzyme in the cholesterol biosynthetic pathway. Alternatively spliced transcript variants encoding different isoforms have been identified for this gene.
------	---

基因ID	23541
基因名	SEC14L2
Swiss	O76054
别名	SEC14L2;C22orf6;SPF;TAP;TAP1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn