

VIPR1 Rabbit pAb

货号: B20072

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	28kDa/46kDa/47kDa/51kDa/55kDa
实测分子量	51kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 31-150 of human VIPR1 (NP_004615.2).
序列	ARLQEECDYVQMIEVQHKQCLEEAQLENETIGCSKMWDNLTCWPATPRGQVVVLACPLIFKLFSSIQGRNVSRSCTDEGW THLEPGPYPIACGLDDKAASLDEQQTMFYGSVKTYTIGY

靶点信息

研究背景	This gene encodes a receptor for vasoactive intestinal peptide, a small neuropeptide. Vasoactive intestinal peptide is involved in smooth muscle relaxation, exocrine and endocrine secretion, and water and ion flux in lung and intestinal epithelia. Its actions are effected through integral membrane receptors associated with a guanine nucleotide binding protein which activates adenylate cyclase. Several transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	7433
基因名	VIPR1
Swiss	P32241
别名	VIPR1;HVR1;II;PACAP-R-2;PACAP-R2;RDC1;V1RG;VAPC1;VIP-R-1;VIPR;VIRG;VPAC1;VPAC1R;VPCAP1R

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn