

KREMEN1 Rabbit pAb

货号: B20073

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	50kDa/51kDa/53kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,22Rv1,A375,Mouse lung,Rat kidney
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 180-380 of human KREMEN1 (NP_001034659.2).
序列	YGEAASTECNSVCFGDHTQPCGGDGRILFDTLVGACGGNYSAMSSVVYSPDFPDYATGRVCYWITRVPGASHIHFSFPLFDIRDSADMVELLDGYTHRVLARFHGRSRPPLSFNVSLDFVILYFFSDRINQAQGFVLYQAVKEELPQERPAVNQTVAEVITEQANLSVSAARSSKVLYVITTPSPHPPQTVPGWTVYG

靶点信息

研究背景	This gene encodes a high-affinity dickkopf homolog 1 (DKK1) transmembrane receptor that functionally cooperates with DKK1 to block wingless (WNT)/beta-catenin signaling. The encoded protein is a component of a membrane complex that modulates canonical WNT signaling through lipoprotein receptor-related protein 6 (LRP6). It contains extracellular kringle, WSC, and CUB domains. Alternatively spliced transcript variants encoding distinct isoforms have been observed for this gene.
------	--

基因ID	83999
基因名	KREMEN1
Swiss	Q96MU8
别名	KREMEN1;ECTD13;KREMEN;KRM1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn